

Explore how our **[NEW cGMP gRNA manufacturing service](#)** can accelerate your **[CRISPR therapeutics project](#)**

EN 

[QUICK ORDER](#)

[SIGN IN](#)



Search f

Search

0

[MY
CART](#)

[Products &
Services](#)

[Explore
Solutions](#)

[Support &
Education](#)

[Tools](#)

[Company](#)

[Contact
Us](#)

Frequently asked questions

Our team has assembled a list of frequently asked questions to help you find answers quickly. Filter using one or more categories to focus on specific topics, or use the search bar to perform a text search.

Search all FAQs:

Can xGen™ UDI-UMI Adapters be used for RNA-seq?

The [xGen™ RNA Library Prep Kit](#) includes a stubby adapter in the workflow so the stocked UDI-UMI adapters are not compatible with those kits.

For customers using homebrew, or other commercial RNA-seq workflows that utilize TA-ligation, these UDI-UMI adapters are compatible. However,

it is important to note that stranded RNA-seq libraries constructed using [xGen UDI-UMI Adapters](#) have the insert strandedness flipped (i.e., to the opposite strand) when compared to standard Illumina TruSeq™ adapters. As such, this strand flip must be handled during data analysis.

Disclaimer: TruSeq™ and Nextera™ are registered trademarks of Illumina, Inc., used with permission. All rights reserved.

Categories:

NEXT GENERATION SEQUENCING

Tags:

NGS PRIMERS XGEN PRIMER NGS ADAPTER ADAPTER

NEXT GENERATION SEQUENCING

Live Chat

